


Solar Absorptance and Total Hemispherical Emittance Measurement System

- Portable Solar Absorptance
Measurement System PM-A2
- Total Hemispherical Emittance
Measurement System PM-E2

Portable Solar Absorptance Measurement System PM-A2

Portable visible spectroscopy with direct driven grating



Measurement Principle

$$\alpha_S(\theta) = \frac{\int_{2500}^{250} \{1 - R(\lambda, \theta)\} I_S(\lambda) d\lambda}{\int_{2500}^{250} I_S(\lambda) d\lambda}$$

R	Spectral reflectance
λ	Wavelength 250~2500nm
θ	Incident angle @7°

Characteristics

- Shortening of the measurement time.
- Measurable from small size sample to large size panel.
- Simple operation, small size, light weight, and cheap.
- High accuracy measurement.
- High-integrity data acquisition system.
- Evaluation data by optical constant method.

Main Specification

System	
Wavelength range	250~2500nm with ozone lamp
Light source	75W Xenon lamp
Measuring methods	<ul style="list-style-type: none"> • Reflection measurement • Transmission measurement • Solar absorptance
Measurement accuracy	<ul style="list-style-type: none"> • Within $\pm 2\%$, • Repeatability: Within $\pm 1\%$ • Solar absorptance ± 0.02
Measuring time	Around 2 minutes/1 scan
Standard sample	<ul style="list-style-type: none"> • Diffuseness evaluation sample (SRS99/Labsphere) • Specular evaluation sample (TFA-25C05-4 Al Deposited mirror/Sigma Koki)
Power source	AC100V $\pm 10\%$ 50/60Hz 250VA
External dimensions	W335 × D375 × H195mm (Excluding outshoot such as integrating sphere)
Weight	Approx. 15kg
Usage environment	<ul style="list-style-type: none"> • Temperature : 15~40°C • Humidity : 10~45%RH
Storage environment	<ul style="list-style-type: none"> • Temperature : 10~45°C • Humidity : 10~50%RH
Accessories & attachments	<ul style="list-style-type: none"> • Instruction manual, • CD for measuring applications and driver installation • USB cable 1m, • Power cable

Spectroscope	
Optical arrangement	Crossed Czerny-Turner arrangement
Focal point distance	134mm
Open area ratio	F/3.4
Diffraction grating	3 kinds included ① 250nm-1200 number/mm ② 800nm-600 numbe/mm ③ 1500nm-300 numbe/mm
Resolution	① < 11.2nm (50-600nm) ② < 22nm (600-1050nm) ③ < 44.5nm (1050-2500nm)
Dispersion	Maximum 15nm/mm
Wavelength drive	Direct driven system
Slit size	φ 2mm (Common in both incoming and outgoing)

Integrating spehere - Data acquisition system with PC	
Internal diameter	φ 100mm
Material	Spectralon®
Detectors used	① Si photodiode × 2 ② InGaAs photodiode × 1 ③ InGaAsPIN photodiode (one step cooling) × 1
Data acquisition system	<ul style="list-style-type: none"> • Data table of spectral reflectance and transmittance • Graphical representation of spectral reflectance and transmittance • Calculation of solar absorptance

Measurement Method

Small size sample

Large size sample




Sample




Uncertainty αs

Standard sample (SRS99)	$\pm 1.0 \sim \pm 2.0\%$
Influence by the multiple reflection of the sample	0.01~0.04





Measurement Example-1

Spectral Transmittance



Measurement Example-2

Spectral Reflectance



Measurement Example-3

Solar Absorptance

Sample	PM-A2	Publication data*
SRS10	0.90	0.89 (0.89)
SRS99	0.01	—
Black Kapton	0.93	0.93 ~ 0.94 (0.94)
Deposited Al Upilex25R	0.39	0.29 ~ 0.33 (0.30)
Deposited Al Upilex25Rnew	0.40	—
CSi100	0.19(Al)	0.17(Al) ~

* JAXA Tukuba, JAXA isas, Kyushu Int.(data)

Delivery Performances

1	JAXA isas	
2	JAXA Mitaka	2011
3	Hokkaido University	2013
4	Nagoya University	2014